ORIGINAL

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TRANSLATION OF THE ORIGINAL USER MANUAL

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EN

The attendant control for the SOLO has been successfully tested according to the standards DIN EN 12182 and DIN EN 12184 at the Technical University of Berlin, Berlin-Cert.

People with visual and reading impairments or intellectual disabilities should print out the user manual available on the website in enlarged letters or have this manual read to them, if necessary. On our homepage, you will also find all safety-related changes which have been made to the user manual. All versions with details of the production period of the product are obtainable at:

www.aat-online.de
- Downloads
- User manuals
- SOLO

The product is state of the art and complies with the associated technical safety requirements. Nevertheless, a low level of residual risks cannot be ruled out. Please observe the safety instructions to this end.

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1 Symbols and safety instructions 04

2 Description of the product 08

3 Mounting the attendant's control 09
3.1 Initial Operation 10
3.2 Driving with the SOLO 12
3.3 Transfer to the control unit 14
3.4 Turn off the SOLO 15

4 Warning notices, malfunction and trouble shooting 15

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EN

1 Symbols and safety instructions

WARNING - Can lead to serious injuries or death!

CAUTION - Can lead to light injuries!

NOTE - Can lead to material damage!

» Before initial operation, read the user's manual carefully and take note of the warning and safety instructions as well as the wheelchair's limiting values.
» Before initial operation, check the electronic function, anti-tilt supports, brackets and, if applicable, the powered wheels (7 bar) and the steering rollers' air pressure.
» Before initial operation, store the SOLO at least one hour in a non-freezing environment.
» Take note of the technical data.
» The approved total weight may not be exceeded.
» Due to the mounting of the SOLO the measurements of the wheelchair can change.
» Operation only by one trained person. This person should be physically and intellectually capable to handle the SOLO safely.
» Practice outside on level ground and at the lowest possible speed.
» For the transportation of persons, not of loads.
» Secure persons who can not seat on their own with a belt.

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- » Never touch rotating or moving parts.
- » If the drive elements block or the driving characteristics change, turn off the SOLO immediately and stop using it.
- » Avoid using the SOLO in heavy rain, ice, snow and on unsecured grounds (sand, mud, loose gravel).
- » Driving on public roads, observe traffic regulations.
- » When crossing streets, intersections and train tracks, use extra caution.
- » Never use the SOLO without anti-tilt supports.
- » Do not drive backward to a wall with mounted anti-tilt supports.
- » Use the lowest possible speed when driving on slopes.
- » Avoid driving down slopes manually with powered wheels.
- » Do not use the SOLO on escalators, stairs or moving belts.
- » Do not drive too closely to stairs or precipices.
- » Do not drive parallel to steep mountain sites or tracks.
- » Keep a safe distance to curbs.
- » Always go backwards over curbs/obstacles.
- » Do not use in connection with hand bikes or other powered wheelchair additions.
- » When not in use, e.g. while using mass transportation (bus, subway, etc.) or when you change seats, turn off the SOLO and engage the locking brakes on the wheelchair.
- » The wheelchair's locking brakes must be adjusted properly.

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EN

- » Turn off the SOLO for transportation.
- » During transportation in vehicles, secure the SOLO according to the relevant legal regulations.
- » Do not use the SOLO plus wheelchair as a seat in a vehicle.
- » Charge the battery pack after each use.
- » Never carry the charger at the charging cable.
- » When charging the SOLO, make sure the cable lies on safe ground.
- » Should the battery pack catch fire, put it out with sand only.
- » Never short-circuit or connect the battery pack's poles.
- » Should the battery pack be damaged or defect, do not use it! Contact the manufacturer or authorized personnel.
- » Keep the plug-in contacts dry.
- » Protect the SOLO against humidity.
- » Protect the charger and the voltage converter for the car from oils, grease, aggressive cleaning agents, thinners or damage.
- » Clean the wheelchair according to the manufacturer's specifications.
- » Do not use a high-pressure cleaning device.
- » We recommend a safety check every two years.
- » Repairs and safety checks may only be done by the manufacturer or authorized personnel.
- » The surface temperature could increase due to the influence of external sources like sunlight.
- » Keep the SOLO away from ignition sources.
- » The SOLO may interfere with other electronic devices' sensitive electromagnetic fields (e.g. anti-theft systems in department stores).

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EN


- Do not operate the SOLO in close proximity to medical devices with life-supporting functions as well as diagnostic devices.
- Only use accessories by AAT Alber Antriebstechnik GmbH.
- Only use the charger of AAT Alber Antriebstechnik GmbH.
- Store the user‘s manual at an appropriate place.

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
3 Mounting the attendant's control

WARNING - Locking bolts must catch!

Step 1:
Remove inserted star shaped grip from the handle adaptor by turning it.

Step 2:
Insert the handle into the handle adaptors and reattach the star shaped grip by tightening.

Repeat steps 1 and 2 on the other side.
> Pushing handle left
> Operating handle right

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EN

Display “wait“
» Connection is established.

Contirm the transfer to the control unit with “+“ or “-“.
» Press “+“ key to permit the transfer to the attendant’s control.
» Press ”-“ to suppress the transfer to the attendant’s control.

Display control unit after transfer.

Standard display control unit
» Battery capacity
» Speed selection

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3.2 Driving with the SOLO

WARNING - Never touch rotating or moving parts!

CAUTION - If the drive elements block or the driving characteristics change, turn off the SOLO immediately and stop using it!

CAUTION - When driving uphill the control unit reacts more sensitively. Please briefly let go of the control grip when the deflective movements are too strong.

NOTE - Due to a longer braking distance the maximum speed forward in combination with the attendant control is 3.6 km/h (driving level 6).

Note the chapter on warnings, error indications and trouble shooting in the original user's manual! The display's light turns itself off after 20 seconds of non-use. During longer brakes of non-use, the automatic switch-off turns the SOLO off after 10 minutes.

Check battery capacity.
» 5 bars = full
» 3 bars = half full
» 1 bar = empty -> charge!

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Move atttandant control in the direction you want to go.

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Press rocker switch in the direction you want to go.

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Panic switch
Press panic switch all the way -> SOLO goes into freewheel. When released, the SOLO brakes immediately.
» This is independent of the travel direction.

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3.3 Transfer to the control unit

Press the ON/OFF key briefly (< 2 seconds).

Confirm the transfer to the control unit with “+“ or “-“.
» Press the “+” to permit the transfer to the control unit!
» Press the “-“ key to suppress the transfer to the control unit!

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4 Warning notices, malfunction and trouble shooting

Please refer to the original manual of the SOLO/SOLO+.

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Translation of the Original EU Declaration of Conformity

We hereby declare under our sole responsibility

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Ehestetter Weg 11

72458 Albstadt/Germany

SRN: DE-MF-000007581

that the wheel hub drives as described in the following for manual wheelchairs

as of serial number: T2001230032
(Basic-UDI-DI) 4054959T2001NK
Risk class: Class I (in accordance with Annex VIII, Rule 1 and 13)

Article no. 335350

Wheel hub drive SOLO+ BLACK EDITION T2001

combined with the optional accessory "Attendant control"

complies with the basic safety and performance requirements and applicable parts of the following regulations and standards.

2017/745/EU	ISO 7176-21; ISO 7176-14; EN 12184; EN 12182; EN ISO 10993-1 IEC 60335-2-29 + A1 & A2 in conjunction with IEC 60335-1
2011/65/EU	EN IEC 63000

The products are marked with CE in compliance with the regulations.

This declaration is only valid for the products with the article numbers indicated above.

This declaration is replaced whenever there is a significant change to the product, product family, standards and laws.

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