

Order form tauron|dts Shower commode chair XXL

Status 2024-01 / Version 2.0



Shaping mobility.

Completing the form on paper

Simply print out the order form, check the relevant component boxes to configure the product, sign the form and send it to us at:

Fax: +49 7248.9186 - 86

Completing the form on a computer or tablet

If you open the document in Acrobat, a menu bar highlighted in colour appears at the top edge.

If the option „Highlight Existing Fields“ is activated, all writeable and selectable fields will be displayed.

- ☐ Circular selection fields indicate that a selection has to be made within the chosen selection.
- ☒ Circular filled selection fields indicate that these options are included in the standard.
- ☐ Square selection fields are not subject to any conditions and are freely selectable and deselectable.

Save the order form locally on your computer under a free selectable name and send it to our customer service via email / fax :

info@dietz-group.com


fax: +49 7248.9186 - 86

Printing errors, mistakes and changes to prices and products reserved. © DIETZ GmbH REHA-PRODUKTE, Karlsbad
This catalogue may be reprinted, in part or in full, only subject to the prior written consent of DIETZ GmbH, Karlsbad.

You can find our general terms and conditions at: <https://shop.dietz-group.de/allgemeine-geschaeftsbedingungen-agb>

Details to your order:

Healthcare supply store		
Customer no.		
Order	Offer based on selling price	Offer based on purchase price
Name of patient/commission		
Delivery address		
Comments		
Company stamp		
Date, signature		



tauron|dts®

	SIDE PANEL STANDARD	SIDE PANEL ABDUCTION WITH STAND-UP AID	SIDE PANEL SWIVELLING
Colours	Aluminium frame lacquered with green colour accents		
Seat width (mm)	560 / 630 / 700		
Seat width front (mm)	-	710 / 780 / 850	-
Seat depth (mm)	450 / 510		
Seat height (mm)	450 – 600		
Overall length (mm)	705 – 747	770 – 806	705 – 747
Overall width (mm)	SW 560: 750 – 792 SW 630: 820 – 862 SW 700: 890 – 932	SW 560: 810 SW 630: 885 SW 700: 960	SW 560: 750 – 792 SW 630: 820 – 862 SW 700: 890 – 932
Overall height (mm)	860 – 1060 (in 25 mm steps)		
Seat cushion (W × D × H)	670 × 450 × 32		
Back cushion (W × D × H)	420 × 70 × 245		
Adjustable back rest angle: (°)	360		
Back height (mm)	410 / 435 / 460		
Side panel height (mm)	200		
Side panel height front (mm)	-	330	-
Armrest length (mm)	580	680	580
Drive-over height (mm)	343 – 493		
Max. load (kg)	250		
Weight (kg)	14	14,5	15


Product details – basic configuration

- Max. load 250 kg
- 3 seat widths adjustable
- 2 seat depths adjustable
- Seat height adjustable in 7 steps of 25 mm
- Backrest height adjustable and 360° angle adjustable
- Upholstered seat and back
- Rubber feet for secure standing
- foldable
- Tested according to the state of the art - > MDR
- Modern look with green colour accents

Price for basic configuration (in EUR)

from 1.390,00


Obesity typology normally distinguishes between A body types and B body types



Typ A

Body type A

- Predominantly male.
- Accumulates fat primarily in the abdominal region, usually beneath the abdominal muscles.
- Normally forces the person to sit with the legs very far apart. It is also usually more comfortable to lean back in order to ease the burden on the abdominal region.




Typ B

Body type B


- Predominantly female.
- Accumulates fat primarily around the bottom and on the upper and lower thighs, usually above the muscles.
- Bodily volume in the bottom region hinders a comfortable sitting posture

Choose the right side panel for every types




Standard side element

These side elements are designed to meet the general requirements of everyday care settings.



Abduction side element with stand-up aid

The abduction side elements add an extra 150 mm onto the seat width so that those who need considerably more space towards the front can sit comfortably on the chair. The side elements curve up and round at the end closest to the user's hands for greater security when getting in and out of the chair.




Swivelling side element

The swivelling side elements on both sides enable easy transfer, e.g. from a wheelchair into the shower commode chair, and vice versa.



Key dimensions and specifications

Together with the height and weight of the user, these are the key dimensions allowing us to align the shower commode chair XXL perfectly to the needs of the user.

Typ A




Typ B



	DATA
	Age of customer
	Height (mm)
	Weight (kg)
1	Seat width back (mm)
2	Seat width front (mm) Only for side panel abduction with stand-up aid
3	Seat depth (mm) Front edge backrest to front edge seat cushion
4	Seat height (mm)

tauron|dts order form

Description



SELECT SIDE PANEL

Basic model		tauron dts commode chair XXL	
	Item no.	720000	
	EUR*	1.390,00	

Side panel	Side panel standard	Side panel abduction with stand-up aid	Side part swivelling	
	Item no.	720100	720101	720102
	EUR*	0,00	280,80	561,60

Selection	Item no.	Description	EUR*
-----------	----------	-------------	------

TOILET BUCKET

720107	Toilet Bucket round (W × H: 300 × 160 mm / Ø 300 mm)	17,40
720108	Toilet Bucket square (L × W × H: 410 × 300 × 160 mm)	53,90

CUSTOM MADE SOLUTIONS

932151	Addition load capacity of max. user weight up to 300 kg	on request
932152	Addition load capacity of max. user weight up to 350 kg	on request

Can only be processed in conjunction with a fully completed customized production report

SERVICE

1015	Service fee for testing and documentation (not discountable)	572,00
------	--	--------

Product group: XXL-Products

Product: tauron|dts

DIETZ GmbH

Reutäckerstraße 12

76307 Karlsbad

info@dietz-group.de

www.dietz-group.de