

Order form CANEO_XL

State 2022-01 / Version 8.0.0



Shaping mobility.

Completing the form on paper

Simply print out the order form, check the relevant component boxes to configure the product, sign the form and send it to us at:

Fax: +49 7248.9186 - 86

Completing the form on a computer or tablet

If you open the document in Acrobat, a menu bar highlighted in colour appears at the top edge. If the option „Highlight Existing Fields“ is activated, all writeable and selectable fields will be displayed.

- Circular selection fields indicate that a selection has to be made within the chosen selection.
- Circular filled selection fields indicate that these options are included in the standard.
- Square selection fields are not subject to any conditions and are freely selectable and deselectable.

Save the order form locally on your computer under a free selectable name and send it to our customer service via email / fax :

info@dietz-group.com

fax: +49 7248.9186 - 86

Printing errors, mistakes and changes to prices and products reserved. © DIETZ Gmbh REHA-PRODUKTE, Karlsbad
This catalogue may be reprinted, in part or in full, only subject to the prior written consent of DIETZ GmbH, Karlsbad.

You can find our general terms and conditions at: <https://shop.dietz-group.de/allgemeine-geschaeftsbedingungen-agb>

CANEO_XL

Folding wheelchair for heavy loads



TECHNICAL	CANEO_XL	
Versions	170 kg	200 kg
Seat width (mm)	560 / 600	
Seat depth* (mm)	400 / 440	
Seat height (mm)	440 / 470	
Back height (mm)	430	
Armrest height	230	
Overall dimensions (W × H × L) (mm)	760/800 × 920-950 × 1150-1200	
Seat angle (degree)	3°	
Back angle (degree)	3° with lumbar curve	
Height of push handles (mm)	920 at seat height 440	
Overall width (mm)	seat width + 200 mm	
Overall width folded (mm)	340	
Overall length without leg rest (mm)	850	
Turning radius min.-max. (mm)	866 - 893	
Armrest height (mm)	230	
Armrest width (mm)	60	
Armrest length (mm)	330 - 380	
Lower leg length (mm)	420 - 510 in steps of 15 mm	
Rear/drive wheels (inch)	24 × 3/8 heavy duty	
Front/swivel wheels (inch)	7 × 1 3/4	
Tipping stability (degree)	10°	
Max. load capacity (kg)	170	200
Weight (kg)*	25	
Materials	Steel frame, high-quality plastic covers	

Product details

Basic configuration

- Tipping stability up to 10°
- Backrest angle 3°, seat angle 3°
- Adjustable backrest (Velcro)
- Reinforced steel sideguards with connection to the sideframe, extended forwards, swiveling to the rear
- Steering wheels shifted more forward and outward, double crossbrace
- Reinforced drive wheel connection
- Lower leg length adjustable in 15 mm steps from 420 up to 510 mm
- Legrests swiveling sideways / removable
- Aluminium foot board
- Park brake/ drum brake (depending on version)
- Drive wheel PU 24" heavy duty , incl. quick release axle, swivel wheel PU 7"
- Reinforced, high-quality plastic covers
- The CANEO_XL is tested according to ISO standard 7176-19 with a power nod system (ISO 75078-2).

*exemplary weight for seat width 60 mm with leg rests, without drum brake

Order form CANEO_XL




Details of your order:

Healthcare supply store		
Customer no.		
Order	Offer based on selling price	Offer based on purchase price
Name of patient/commission		
Delivery address		
Comments		
Company stamp		
Date, signature		

Order form CANEO_XL

Key dimensions and specifications

Together with the height and weight of the user, these are the key dimensions allowing us to align the wheelchair perfectly to the needs of the user.



	DATA
	Age of customer
	Height (cm)
	Weight (kg)
1	Seat width (cm)
2	Seat depth (cm) Length of upper leg less two finger breadths
3	Height of upholstered backrest (cm) From seat surface to below the shoulder blades
4	Lower leg length (cm) For the correct height of the foot plate
5	Seat height, front (cm) Length of lower leg + sufficient clearance between foot plate and ground (at least 4 cm)
6	Seat height, back (cm) Depends on front seat height, seat angle and seat depth

Order form CANEO_XL

Designation

SELECT THE RIGHT SEAT WIDTH

	170 kg	200 kg	170 kg	200 kg
Seat width (cm)	56	56	60	60
With locking brake				
Item no.	151161	151056	151162	151060
price*	2.125,00	2.525,00	2.225,00	2.625,00
With drum + locking brake				
Item no.	151163	151156	151164	151160
price*	2.325,00	2.725,00	2.425,00	2.725,00

No tools are needed to reconstruct the seat depth from 40 cm to 44 cm - independent of seat with.


Selection	Item no.	Description	EUROS*
-----------	----------	-------------	--------

SELECT THE SEAT HEIGHT

Standard	Seat height 47 cm	assembled	0,00
999944	Seat height 44 cm	assembled	25,00

The modification can be done by your specialist retailer without additional material.

SELECT LEGREST




Standard	Standard legrest ¹ / left, foot plates height-adjustable	assembled	0,00
Standard	Standard legrest ¹ / right, foot plates height-adjustable	assembled	0,00
198004	Amputee legrest, left	assembled	121,00
198000	Amputee legrest, right	assembled	121,00
198010	Elevating legrest / left, angle adjustable (SW 56 – 60 cm)	assembled	179,00
198016	Elevating legrest / right, angle adjustable (SW 56 – 60 cm)	assembled	179,00

¹ The standard legrest is delivered with aluminium footboard.

Order form CANEO_XL

Selection	Item no.	Description		EUROS*
-----------	----------	-------------	--	--------

OPTIONS OF POSITIONING



Standard	Adjustable back (Velcro)	assembled	0,00
----------	--------------------------	-----------	------

Designation	Seat width
-------------	------------

ANTI-DECUBITUS

Seat width (cm)	56			60		
	Selection	Item no.	Price	Selection	Item no.	Price
Vicair XXtra O2* cushion, black Seat depth 45 cm height 10 cm		420220	950,00		420221	950,00

*Vicair XXtra O2 Anti-Decubitus cushion is breathable. The padding consists of Smart Cells which minimizes slide- and shear forces. The measurements are tested on the CANEO_XL - variations towards the measurements of the wheelchair of +/-1 cm are negligible.

WHEELCHAIR SET



Walking stick holder 2 sticks



Lap belt



Anti-tip wheels

404076	Wheelchair set consisting of: Walking stick holder for 1-2 sticks, anti-tip wheels right and left, lap belt with 2-point seat belt lock „Press“		133,60
--------	---	--	--------

Order form CANEO_XL

Selection	Item no.	Description	Comment	EUROS*
-----------	----------	-------------	---------	--------

BASIC MODEL OPTIONS ACCORDING TO SEAT WIDTH

196252	Headrest incl. bracket (SW 56 cm)		not suitable for car use	375,00
195045	Activity tray (SW 56 cm)			251,00
195047	Activity tray (SW 60 cm)			284,00
190054	Activity tray half size, foldable down, left / right			196,00

VARIOUS OPTIONS FOR ALL SEAT WIDTH

195017	Power node system with long lap belt & headrest, only possible with Standard seat height 47 cm		Crash test approved ISO 7176-19	1.262,00
196236	Anti-tip - wheel, left		assembled	43,00
196237	Anti-tip - wheel, right		assembled	43,00
190033	Park brake one-hand drive, left		assembled	258,00
190034	Park brake one-hand drive, right		assembled	258,00
700106	Wheel, set air, swivel wheel 7 × 1 3/4"		assembled	37,00
700103	Wheel, set air, drive wheel 24 × 1 3/8" CST		assembled	37,00
198065	Brake lever extension, park brake, left		assembled	24,00
198064	Brake lever extension, park brake, right		assembled	24,00
198069	Spoke guard with fittings, transparent for drive wheel 24"		assembled	88,00
196255	Walking stick holder with bracket for 25 mm cane			18,60
402114	Lap belt (press with metal buckle) 1800 mm adjustable in length		not suitable for car use	29,00
190053	Infusion holder, stainless steel with two arm holder of plastic, ca. 100 cm			152,00
196256	Push rim cover silicon, 24" wheel, black/grooved			64,00
196257	Push rim cover silicon, 24" wheel, black/smooth			64,00
196259	Folding lock black			59,00

Manual wheelchairs

DIETZ Power Electric wheelchairs

Rollators

Walking aids

Scooters

Bath & WC

XXL-Products

DIETZ GmbH
Reutäckerstraße 12
76307 Karlsbad
info@dietz-group.de
www.dietz-group.de